

TEST TUBES (RIMLESS)

Item no.	Size	Product Code	Glass Spec	Qty per Tray
1	10x75	KASD001	9.75 x 0.75	1500
2	12 x 75	KASD002	12.0 x 0.8	1000
3	12 x 100	KASD003	12.0 x 0.8	1000
4	15 x 100	KASD004	15.2 x 0.8	621
5	16 x 100	KASD005	16.0 x 0.9	561
6	16 x 125	KASD006	16.0 x 0.9	561
7	16 x 150	KASD007	16.0 x 0.9	561
8	18 x 150	KASD008	18.5 x 0.95	418
9	20 x 150	KASD009	20.7 x 0.95	320
10	24 x 125	KASD010	24.5 x 1.05	231
11	24 x 150	KASD011	24.5 x 1.05	231
12	18 x 125	KASD012	18.5 x 0.95	418

TEST TUBES (WITH RIM)

Item no.	Size	Product Code	Glass Spec	Qty per Tray
1	16 x 125	KASD013	16.0 x 0.9	1000
2	16 x 150	KASD014	16.0 x 0.9	1000
3	20 x 150	KASD015	20.7 x 0.95	1000
4	24 x 150	KASD016	24.5 x 1.05	1000

TABLET TUBES - POLY TOPS

Item no.	Size	Product Code	Glass Spec	Length	QTY per Tray
1	No.1	TAB - KASD017	15.2 x 0.8	49/50mm long	572
2	No.2	TAB - KASD018	18.5 x 0.95	49/50mm long	396
3	No.3	TAB - KASD019	18.5 x 0.95	59/60mm long	396
4	No.4	TAB - KASD020	24.5 x 1.05	49/50mm long	217
5	No.5	TAB - KASD021	24.5 x 1.05	59/60mm long	217
6	No.6	TAB - KASD022	24.5 x 1.05	74/75mm long	217
7	No.7	TAB - KASD023	24.5 x 1.05	89/90mm long	217