



## Kayla Africa Suppliers & Distributors Cc



Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p><b>Features:</b> Fully stainless steel structure Automatic sterilization,easy to operate Rapid Sterilization in 4-6 minutes European CLASS N</p> <p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134℃ Adjustment of temperature: 105~135℃ Timer: 0~60min Power: 1.5KW / AC220V 50HZ</p>	TS-B24	24L Ø270*H410mm 3pcs sterilizing plate	560x350x500	18
	<p><b>Features:</b> Fully stainless steel structure Automaticsterilization ControlledbyMicro-computer Digital display of working status Rapid Sterilization in 4-6 minutes European CLASS N</p> <p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134℃ Adjustment of temperature: 105~135℃ Timer: 0~99min Power: 1.5KW / AC220V 50HZ</p>	TS-BD24	24L Ø270*H410mm 3pcs sterilizing plate	560x350x500	18

Please find your preferred models and send email to us for request Latest Best Price



## Kayla Africa Suppliers & Distributors Cc



Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog -

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p><b>Features:</b> Sterilizing chamber made of stainless steel Rapid Sterilization in 4-6 minutes Steam-water inner circulation Safety Protection Devices Class N European Standard</p> <p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~60min Power: 1.5KW / AC220V 50HZ / AC 110V 60Hz</p>	TS-AJ20	20L Ø250*H420mm 3pcs sterilizing plate	480x480x384	40
		TS-AJ24	24L Ø250*H520mm 3pcs sterilizing plate	580x480x384	45
	<p><b>Features:</b> Fully Stainless Steel Structure Rapid Sterilization in 4-6 minutes Steam-water inner circulation Safety Protection Devices Class N European Standard</p> <p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~60min Power: 1.5KW / AC220V 50HZ / AC 110V 60Hz</p>	TS-A20	20L Ø250*H420mm 3pcs sterilizing plate	480x480x384	38
		TS-A24	24L Ø250*H520mm 3pcs sterilizing plate	580x480x384	44

Please find your preferred models and send email to us for request Latest Best Price



## Kayla Africa Suppliers & Distributors Cc



Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p><b>Features:</b> Fully Stainless Steel Structure Rapid Sterilization in 4-6 minutes Digital Display of Working status, touch type key Steam-water inner circulation Safety Protection Devices Drying function Class N European Standard</p> <p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~99min Power: 1.5KW / AC220V 50HZ</p>	TS-AD20	20L Ø250*H420mm 3pcs sterilizing plate	480x480x380	40
		TS-AD24	24L Ø250*H520mm 3pcs sterilizing plate	580x480x384	45
	<p><b>Features:</b> Sterilizing chamber made of stainless steel Rapid Sterilization in 4-6 minutes Steam-water inner circulation Safety Protection Devices Class N European Standard</p> <p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~60min Power: 1.5-2.7KW / AC220V 50HZ</p>	TS-C20	20L Ø250*H420mm 3pcs sterilizing plate	480x480x384	38
		TS-C24	24L Ø250*H520mm 3pcs sterilizing plate	580x480x384	44
		TS-C35	35L Ø300*H500mm 4pcs sterilizing plate	730x550x510	73
		TS-C50	50L Ø340*H550mm 4pcs sterilizing plate	770x600x540	86

Please find your preferred models and send email to us for request Latest Best Price



## Kayla Africa Suppliers & Distributors Cc




Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p>Features: Sterilizing chamber made of stainless steel Rapid Sterilization in 4-6 minutes with 3 fixed cycles of sterilizing programs for option or preset by yourself Steam-water inner circulation Whole process is automatic, including add water Safety Protection Devices With drying function Class N European Standard</p> <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~99min Power: 1.5- 2.7KW / AC220V 50HZ</p>	TS-CD35	35L Ø300*H500mm 4pcs sterilizing plate	730x550x510	73
		TS-CD50	50L Ø340*H550mm 4pcs sterilizing plate	770x600x540	86
	<p><b>Features:</b> Pulsating Vacuum Desktop Steam sterilizer. Automatic Sterilization. LCD display working status, touch keys. Vacuum date can reach -0.8bar, fits various instruments, including package, unpackaged, solid, hollow, multi-aperures and inset pipes. Built-in mini printer. Safety Protection Devices . BOWIE&amp;DICK testing is available. Can inactivate HIV, HBV, mad cow virus and backllus. Class B European Standard.</p>	TS-ADDV18	18L Ø250*H360mm 3pcs sterilizing plate	560x470x400	47
		TS-ADDV24	24L Ø250*H470mm 3pcs sterilizing plate	675x470x400	53
	<p><b>Technical Data:</b> Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C; 121/134°C for 18/24L Timer: 0~99min Power: 2-2.7KW / AC220V 50HZ</p>	TS-DDV35	35L Ø300*H500mm 4pcs sterilizing plate	730x550x510	93
		TS-DDV50	50L Ø340*H550mm 3pcs sterilizing plate	770x600x540	106

Please find your preferred models and send email to us for request Latest Best Price



## Kayla Africa Suppliers & Distributors Cc




Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>1. Adopting European N standard, with thermal vacuum and drying functions, it can be used for A-type hollow instruments and solid instruments. Thermodynamic vacuum guarantees elimination of cold air from the cavity and makes sure that reaches the optimal effect of sterilization.</li> <li>2. Specially design for dental clinic field, it could sterilize 12 pieces hand piece at same time, simply and efficiently.</li> <li>3. Sliding door design for easy secure operating-for TS-DDN2</li> <li>4. Microprocessor control, LCD or LED display and modularized panel for easy operation.</li> <li>5. Automatic malfunction detection system makes easy maintaining.</li> </ol> <p><b>Technical Data:</b>            Working pressure: Min. -0.9bar, Max. 2.3 bar            Working temperature: 121 °C, 134 °C            Water consumption per one cycle: 0.12L~0.16L(TS-DDN2); 0.16-0.18L(TS-DDN7/10/12/18/24)            Timer: 0-90 min            Power: 1.2KW-1.7KW / AC220V 50HZ            Note: The price is based on power supply (220V, 50Hz). If need to change power, will cost extra fees</p>	TS-DDN2	2L Ø88*H175mm	350*284*376	26
		TS-DDN7	7L Ø170*H320mm	540*370*360	37
		TS-DDN10	10L Ø192*H320mm	540*370*360	40
		TS-DDN12	12L Ø192*H340mm	560*450*400	43
		TS-DDN18	18L Ø245*H360mm	560*450*400	47
	TS-DDN24	24L Ø245*H470mm	670*450*400	53	

Please find your preferred models and send email to us for request Latest Best Price



## Kayla Africa Suppliers & Distributors Cc








Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>In conformity with EN13060 standard with opening water tank on the top which makes wash and cleaning more convenient</li> <li>It adopts European B standard with 3-times pre-vacuum. Residual humidity is less than 0.2% after vacuum drying.</li> <li>Vacuum rate reaches -0.8bar. Suitable for various instruments, including package, unpackaged, solid, hollow, multi-apertures and inset pipes</li> <li>Fuzzily Micro-computer control, bright digital display. The interface is modularized panel for easy operation</li> <li>With Bowie &amp; Dick test verifying steam penetration into porous bodies; with Helix Test verifying steam penetration into hollow bodies.</li> <li>Optional configurations: Mini printer can be installed to record the process of sterilization. USB flash disk for electronic storage of sterilization data (It is standard configuration for TS-DDBN18/24.)</li> </ol> <p><b>Technical Data:</b> Working pressure: Min. -0.9bar, Max. 2.3 bar Working temperature: 121℃, 134℃ Water consumption per one cycle: 0.12L~0.16L (TS-DDB2); 0.16L~0.18L (TS-DDB7/10/12/16/18/24) Timer: 0-90 min Power: 1.2KW-1.4KW / AC220V 50HZ Note: The price is based on power supply (220V, 50Hz). If need to change power, will cost extra fees</p>	TS-DDB2	2L Ø88*H175mm	350*284*376	26
		TS-DDB7	7L Ø170*H320mm	540x370x360	37
		TS-DDB10	10L Ø192*H320mm	540*370*360	40
		TS-DDB16	16L Ø240*H360mm	445*550*395	47
		TS-DDB12	12L Ø192*H340mm	560*450*400	43
		TS-DDBN18	18L Ø245*H360mm	560*450*400	47
		TS-DDBN24	24L Ø245*H470mm	670*450*400	53

Please find your preferred models and send email to us for request Latest Best Price



## Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston,

Johannesburg Email: [k.a.s.d@telkomasa.net](mailto:k.a.s.d@telkomasa.net)

Website ID: [www.kasd.co.za](http://www.kasd.co.za)

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

### Table Top Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Gross Weight (kg)
	<p><b>Features :</b></p> <ol style="list-style-type: none"> <li>1.LCD display running status</li> <li>2. Multiple security protective device</li> <li>3. Automatic door can effectively prevent scald</li> <li>4. Ultimate vacuum value can reach -90KPa or above</li> <li>5. With build-in steam generator</li> <li>6. Multiple programs include: packed items, unpacked items, custom program, rapid program, BD testing program, vacuum testing program, preheat program and drying program.</li> </ol> <p><b>Technical Data:</b></p> <p>Design pressure: -0.1/0.3Mpa            Working temperature: 105-136 °C            Timer: 0-99mins            Power: 2kw-7kw/AC220V/50Hz</p>	TS-DDZ18	18L Ø250*H350mm	585*515*450	55
		TS-DDZ24	24L Ø250*H450mm	685*515*450	65
		TS-DDZ45	45L Ø316*H618mm	830*640*550	115
		TS-DDZ60	60L Ø385*H595mm	1002*730*560	160
		TS-DDZ80	80L Ø385*H717mm	1002*730*560	180

Please find your preferred models and send email to us for request Latest Best Price