




Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston, Johannesburg

Email: k.a.s.d@telkomasa.net Website ID: www.kasd.co.za

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Automatic control in processes of filling water, heating, sterilization, steam e*haust and dryness. 3. Door safety lock system 4. Built-in water tank 5. Preheating controlfunction 6. Steam-water inner circulation system 7. Protection of water lacking, over pressure, over heating <p>Technical Data:</p> <p>Working pressure: 0.22Mpa Working temperature: 134 °C Adjustment of temperature: 105~134 °C Timer: 0~99min Power:2.5-4.5KW / AC 220V±22V, 50Hz±1Hz</p>	VA-SA35	35L Ø318*H450mm 2pcs sterilizing baskets	600*410*1030	77
		VA-SA50	50L Ø340*H550mm 2pcs sterilizing baskets	650*410*1140	83
		VA-SA75	75L Ø400*H600mm 2pcs sterilizing basket	690*470*1140	98
		VA-SA100	100L Ø440*H650mm 2pcs sterilizing baskets	730*510*1270	110

Please find your preferred models and send email to us for request Latest Best Price



Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston, Johannesburg

Email: k.a.s.d@telkomasa.net Website ID: www.kasd.co.za

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel material 2. Automatic e*hausting cool air, heating and sterilization 4 foot wheels for easy movement Hand wheel type of quick-open door <p>Double scale indication pressure gauge</p> <ol style="list-style-type: none"> 5. With PP/TT test connector <p>Technical Data:</p> <p>Working pressure: 0.22Mpa Working temperature: 134 °C Adjustment of temperature: 105~134 °C Timer: 0~60min Power: 2.5-4.5KW / AC220V 50HZ</p>	VA-SJ35	35L Ø318*H450mm 2pcs sterilizing baskets	450*450*1010	56
		VA-SJ50	50L Ø340*H550mm 2pcs sterilizing baskets	510*470*1130	68
		VA-SJ75	75L Ø400*H600mm 2pcs sterilizing baskets	560*560*1120	80
		VA-SJ100	100L Ø440*H650mm 2pcs sterilizing baskets	540*560*1250	85

Please find your preferred models and send email to us for request Latest Best Price




Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston, Johannesburg

Email: k.a.s.d@telkomasa.net Website ID: www.kasd.co.za

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features: 28 Liter VerticalType Autoclave Full automatic control by computer With touch keys, digital display Hand wheel type quick open door Suitable for fast 4~6 min sterilization Full stainless steel material Over temperature and pressure protection: when inside chamber pressure over 0.24 MPa Water lacking protection Automatically shut off with beep after sterilization finished</p> <p>Optional: Drying Function (USD85.00 if need)</p> <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134 °C Adjustment of temperature: 105~134 °C Timer: 0~99min Power: 2KW / AC220V 50HZ</p>	<p>VA-SD28</p>	<p>28L Ø270*H500mm 2pcs sterilizing baskets</p>	<p>400*370*700</p>	<p>23</p>

Please find your preferred models and send email to us for request Latest Best Price




Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston, Johannesburg

Email: k.a.s.d@telkomasa.net Website ID: www.kasd.co.za

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features: Vertical Pressure Steam Autoclave Digital display, touch key Hand wheel type of quick-open door Automatic e*hausting cold air, heating and sterilization 4 foot wheels for easy movement Double scale indication pressure gauge Safety Protection Devices</p> <p>Optional: Drying Function (USD85.00 if need)</p> <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134 °C Adjustment of temperature: 105~134 °C Timer: 0~99min Power: 2.5-4.5KW / AC220V 50HZ</p>	VA-SD35	35L Ø318*H450mm 2pcs sterilizing baskets	450*450*1010	56
		VA-SD50	50L Ø340*H550mm 2pcs sterilizing baskets	510*470*1130	68
		VA-SD75	75L Ø400*H600mm 2pcs sterilizing baskets	560*560*1120	80
		VA-SD100	100L Ø440*H650mm 2pcs sterilizing baskets	540*560*1250	85

Please find your preferred models and send email to us for request Latest Best Price



Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston, Johannesburg

Email: k.a.s.d@telkomasa.net Website ID: www.kasd.co.za

Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Micro-computer control, LCD display of working status, touch key 3. Fully automatic route for water inlet, pulse vacuum, heating, sterilizing, draining, e*hausting, vacuum drying. 4. Three times pulse vacuum 5. With 4 fi*ed programs for applicance, medical dressing, rubbers and fluid. 6. PT/TT testing connector 7. Equipped with 0.22 µm bacteria filter 8. Steam-water inner circulation system <p>Optional: Printer function (USD335.00 if need)</p> <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134 °C Adjustment of temperature: 0~134 °C Timer: 1~99min Power: 1.5-4.5KW / AC220V 50HZ</p>	VA-SV35	35L Ø318*H450mm 2pcs sterilizing baskets	650*410*1140	95
		VA-SV50	50L Ø340*H550mm 2pcs sterilizing baskets	690*470*1140	105
		VA-SV75	75L Ø400*H600mm 2pcs sterilizing baskets	730*510*1270	125
		VA-SV100	100L Ø440*H650mm 2pcs sterilizing baskets	770*550*1350	140

Please find your preferred models and send email to us for request Latest Best Price



Kayla Africa Suppliers & Distributors Cc

Address: Unit C9, Wadeville Business Park, Cnr Steenbras & Shad Roads, Wadeville, Germiston, Johannesburg

Email: k.a.s.d@telkomasa.net Website ID: www.kasd.co.za


Tel: +2711 027 7221, Mobile / WhatsApp : +27 83 488 5196

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel 2. Over temperature & over pressure auto-protection 3. Double scale indication pressure gauge <p>Safe protection of water lacking</p> <ol style="list-style-type: none"> 4. Automatically shut off with beep reminding after sterilization <p>Technical Data:</p> <p>Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~134 °C Timer: 0~60min Power: 2.5-6KW / AC220V 50HZ</p>	VA-FJ35	35L Ø318*H450mm 2pcs sterilizing baskets	480*460*850	42
		VA-FJ50	50L Ø340*H550mm 2pcs sterilizing baskets	520*520*980	50
		VA-FJ75	75L Ø400*H600mm 2pcs sterilizing baskets	560*560*980	70
		VA-FJ100	100L Ø440*H650mm 2pcs sterilizing baskets	590*590*1080	85
		VA-FJ120	120L Ø480*H660mm 2pcs sterilizing baskets	600*640*1140	100
		VA-FJ150	150L Ø510*H740mm 2pcs sterilizing baskets	670*690*1140	110

Please find your preferred models and send email to us for request Latest Best Price

Vertical Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full Stainless steel 2. Digital display of working status, touch key 3. Discharge cool air automatically, and steam discharging automatically after sterilization. 4. Drying function optional by customer's request <p>Optional: Drying Function (USD85.00 if need)</p> <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~134°C Timer: 0~99min Power: 2.5-6KW / AC220V 50HZ</p>	VA-FD35	35L Ø318*H450mm	480*460*850	42
		VA-FD50	50L Ø340*H550mm	520*520*980	50
		VA-FD75	75L Ø400*H600mm	560*560*980	70
		VA-FD100	100L Ø440*H650mm	590*590*1080	85
		VA-FD120	120L Ø480*H660mm	600*640*1140	100
		VA-FD150	150L Ø510*H740mm 2pcs sterilizing	670*690*1140	110

Please find your preferred models and send email to us for request Latest Best Price