

KASD® Sterile Syringe Filter

KASD® Sterile syringe filters are available with Polyethersulfone (PES), Polyesteramide(Nylon), Cellulose Acetate(CA), Polytetrafluoroethylene(PTFE), Polyvinylidene Fluoride(PVDF). And each filter is individually packed and sterilized by Gamma Radiation. Every Syringe Filter is printed with expiry date for easy QC tracking.

Application

- Tissue culture media preparation
- Sterile filtration and clarification of biological fluids
- Probe solutions
- Protein and enzyme filtration
- Hybridization buffers
- Other aqueous solutions





Kayla Africa Suppliers & Distributors

Product Information

Item	Description	Packaging (pcs/pk)
PES Syringe Filter		
KASD13PES022S	Pore:0.22µm Diameter: 13mm	100
KASD13PES045S	Pore:0.45µm Diameter: 13mm	100
KASD25PES022S	Pore:0.22µm Diameter: 25mm	50
KASD25PES045S	Pore:0.45µm Diameter: 25mm	50
KASD30PES022S	Pore:0.22µm Diameter: 30mm	50
KASD30PES045S	Pore:0.45µm Diameter: 30mm	50
CA Syringe Filter		
KASD13CA022S	Pore:0.22µm Diameter: 13mm	100
KASD13CA045S	Pore:0.45µm Diameter: 13mm	100
KASD25CA022S	Pore:0.22µm Diameter: 25mm	50
KASD25CA045S	Pore:0.45µm Diameter: 25mm	50
KASD30CA022S	Pore:0.22µm Diameter: 30mm	50
KASD30CA045S	Pore:0.45µm Diameter: 30mm	50
PVDF Syringe Filter(Hydrophobic)		
KASD13PVDF022HS	Pore:0.22µm Diameter: 13mm	100
KASD13PVDF045HS	Pore:0.45µm Diameter: 13mm	100
KASD25PVDF022HS	Pore:0.22µm Diameter: 25mm	50
KASD25PVDF045HS	Pore:0.45µm Diameter: 25mm	50
KASD30PVDF022HS	Pore:0.22µm Diameter: 30mm	50
KASD30PVDF045HS	Pore:0.45µm Diameter: 30mm	50
PVDF Syringe Filter(Hydrophilic)		
KASD13PVDF022LS	Pore:0.22µm Diameter: 13mm	100
KASD13PVDF045LS	Pore:0.45µm Diameter: 13mm	100
KASD25PVDF022LS	Pore:0.22µm Diameter: 25mm	50
KASD25PVDF045LS	Pore:0.45µm Diameter: 25mm	50
KASD30PVDF022LS	Pore:0.22µm Diameter: 30mm	50
KASD30PVDF045LS	Pore:0.45µm Diameter: 30mm	50
GF Syringe Filter		
KASD13GF070S	Pore:0.7µm Diameter: 13mm	100
KASD13GF120S	Pore:1.2µm Diameter: 13mm	100
KASD25GF070S	Pore:0.7µm Diameter: 25mm	50



Kayla Africa Suppliers & Distributors

KASD25GF120S	Pore:1.2µm	Diameter: 25mm	50
KASD30GF070S	Pore:0.7µm	Diameter: 30mm	50
KASD30GF120S	Pore:1.2µm	Diameter: 30mm	50
Nylon66 Syringe Filter			
KASD13NY022S	Pore:0.22µm	Diameter: 13mm	100
KASD13NY045S	Pore:0.45µm	Diameter: 13mm	100
KASD25NY022S	Pore:0.22µm	Diameter: 25mm	50
KASD25NY045S	Pore:0.45µm	Diameter: 25mm	50
KASD30NY022S	Pore:0.22µm	Diameter: 30mm	50
KASD30NY045S	Pore:0.45µm	Diameter: 30mm	50
PTFE Syringe Filter (Hydrophobic)			
KASD13PTFE022HS	Pore:0.22µm	Diameter: 13mm	100
KASD13PTFE045HS	Pore:0.45µm	Diameter: 13mm	100
KASD25PTFE022HS	Pore:0.22µm	Diameter: 25mm	50
KASD25PTFE045HS	Pore:0.45µm	Diameter: 25mm	50
KASD30PTFE022HS	Pore:0.22µm	Diameter: 30mm	50
KASD30PTFE045HS	Pore:0.45µm	Diameter: 30mm	50
PTFE Syringe Filter (Hydrophilic)			
KASD13PTFE022LS	Pore:0.22µm	Diameter: 13mm	100
KASD13PTFE045LS	Pore:0.45µm	Diameter: 13mm	100
KASD25PTFE022LS	Pore:0.22µm	Diameter: 25mm	50
KASD25PTFE045LS	Pore:0.45µm	Diameter: 25mm	50
KASD30PTFE022LS	Pore:0.22µm	Diameter: 30mm	50
KASD30PTFE045LS	Pore:0.45µm	Diameter: 30mm	50

Note:

Various membrane medium with pp or glass fiber pre-filter are available.