

KASD Moisture Analyser



Technical parameters :

Model : KASD-100MW

Range(g) : 110 110 110

Readability(g) : 0.005 0.002 0.001

Operating Temp.Range : 5°C---35°C

Display : LCD

Calibration : External

Pan size (mm) : Φ 90

Dimensions(W*L*H) : 200 x 180 x 300

Weight(Kg) : 5.8

Heat source : Halogen lamp

Temp.sensor : PT-100

Temp.range : 40---199°C

Temp.step : 1°C

Moisture range : 0.00%---100%

Moisture readability : 0.01%

Dry residual range : 100.00%---0.00%

Dry residual readability : 0.01%

ATRO : 100%---999%

ATRO : 0%---999%

Temperature setting : 40---199°C

Time setting : 1---99min by 1s step

Autostop : 0.1%---9.9% within 1 min

Number of storage : 15

Heater : 220V+15% 50Hz

115V+15%

400W

R232

Balance : External power supply

Output : 9V

Power supply : 220V+15% 50Hz

Selectable parts : Printer, Aluminum case