

Poulten & Graf

SUPERIOR LABORATORY PRODUCTS



Bottle Top Dispenser FORTUNA® OPTIFIX®,
comes with 2 PP adaptors 40+45, 1 PTFE-intake
tube, single pack

OPTIFIX® BASIC

OPTIFIX® SOLVENT

Capacity	grad.
0,4- 2 ml	0,1 ml
1 -5 ml	0,2 ml
2-10 ml	0,2 ml
6- 30 ml	1,0 ml
10 - 50 ml	1,0 ml
20- 100 ml	2,0 ml
40 - 200 ml	5,0 ml
60- 300 ml	5,0 ml
100- 500 ml	5,0 ml

Art.-No.

101.080-27
101.080-33
101.080-38
101.080-44
101.080-48
101.080-51
101.080-53
101.080-56
101.080-58

Art.-No.

101.081-27
101.081-33
101.081-38
101.081-44
101.081-48
101.081-51
101.081-53
101.081-56
101.081-58



Bottle Top Dispenser FORTUNA® OPTIFIX®,
comes with 2 PTFE adaptors 40+45,
1 PTFE-intake tube, single pack

FORTUNA UNIVERSAL® OPTIFIX® HF

OPTIFIX® SAFETY

OPTIFIX® SAFETY S

Capacity	grad.
0,4 - 2 ml	0,1 ml
1 - 5 ml	0,2 ml
2 -10 ml	0,2 ml
6- 30 ml	1,0 ml
10 - 50 ml	1,0 ml
20- 100 ml	2,0 ml

Art.-No.
-
101.106-33
101.106-38
101.106-44
101.106-48
-

Art.-No.
101.107-27
101.107-33
101.107-38
101.107-44
101.107-48
101.107-51

Art.-No.
101.107-27S
101.107-33S
101.107-38S
101.107-44S
101.107-48S
101.107-51S



Made
in
Germany

Only at Poulten & Graf:
Bottle Top Dispenser from 2 - 500 ml!

Poulten & Graf GmbH

Karl-Carstens-Str. 10,
97877 Wertheim, Germany
Tel.: +49 9342 9229-0
Fax: +49 9342 9229-80
www.poulten-graf.de
FORTUNA@poulten-graf.com
www.shop.poulten-graf.de

- OPTIFIX® BASIC: for dispensing of the common chemicals
- OPTIFIX® SOLVENT: for dispensing of solvents
- FORTUNA UNIVERSAL® OPTIFIX® HF: for hydrofluoric acid
- OPTIFIX® SAFETY/SAFETY S: for fuming liquids

When dispensing liquids, do not compromise - OPTIFIX® bottle top dispensers are ALWAYS the best choice!

Our aim is for all your dispensing to be both safe and precise! Therefore we have developed the perfect OPTIFIX® Dispenser for all liquid chemicals.

Five OPTIFIX® types are available. User friendly design and tested precision are further decisive advantages of the OPTIFIX® - dispenser. Virtually indestructible through the use of high quality material the OPTIFIX® dispenser is one of the highlights of our product programm.

The five FORTUNA® OPTIFIX® - models:



OPTIFIX® BASIC

for dosing of the most common chemicals, acids and crystallizing liquids (no HF), with dispensing cylinder made of glass and with PTFE-encased piston.



OPTIFIX® SOLVENT

for dispensing of pure organic solvents (no HF), with dispensing cylinder and piston made of glass, also available.



FORTUNA UNIVERSAL® OPTIFIX® HF

for all liquids that can be handled with the OPTIFIX® dispensers and for hydrofluoric acid (HF), with dispensing cylinder made of PTFE and with solid PTFE-encased piston



OPTIFIX® SAFETY

for dispensing of strong acids and aggressive chemicals (no HF), with dispensing cylinder made of glass and PTFE-encased piston with sealing lip



OPTIFIX® SAFETY S

for dispensing of strong fuming acids and very abrasive liquids (no HF), with safety stopcock and additional safety valve in ventilation tube, with dispensing cylinder made of glass and PTFE-encased piston with sealing lip



Volume Setting System

- fast and easy volume setting
- additional fine adjustment for right- and left hand use
- easy to read scale

Dispensing Piston

- proven technology
- made of glass, solid PTFE-encased, exception: SOLVENT: not PTFE-encased
- virtually indestructible

Valve Block

- high technology optimal function
- 100 % tightness
- no metal springs - no corrosion

OPTIFIX® BASIC and SOLVENT

can be sterilized at
+121°C (+248°F)/
1.04 bar.

Finger Protection Cap

- ergonomic type
- timeless design
- with colour-code-system

Precision Dispensing Cylinder

- made of glass, dosing procedure visible (except type HF, made of PTFE)
- coated with plastic - protection against breakage
- graduation clear readable

FIX Adaptor

- to align the dispenser on a bottle
- made of Polypropylene for BASIC and SOLVENT
- made of PTFE for HF and SAFETY